

A COMPLETE INSTRUMENT TRAY SOLUTION

"The Complete System!"



CBM was founded in 1945, and moved to the current site in 1962.

CBM has now developed over an area of 30,000 m², 10,000 m² of which are under cover.

- High quality designers and engineers who pride themselves on the development and compliance to international standards.
- They are reputed manufacturers of high quality stainless steel baskets and rigid containers.
- At Richards Medical, we have developed a partnership in Australia and New Zealand with CBM that brings the CBM products and experience to Australasian customers.
- Richards Medical is delighted to present the CBM products to you. Experience for yourself the difference of how excellence in design, manufacturing, engineering and sheer attention to detail will benefit your department.
- Don't just talk to us about the products, ask us how we can help you with a replacement program.



MINI BASKETS SUITABLE FOR FINE INSTRUMENTS, INCLUDING PLASTICS, ENT, OPHTHALMICS & DENTAL

"Give your delicate instruments Maximum Protection!"

Small Instrument Basket 170 x 135 x 35 mm

Basket has a Silicone mat on the inner base to secure instruments

Length Approx: 170 mm
Width Approx: 135 mm
Height Approx: 35 mm
Product Code: CBM67S



Small Instrument Basket 270 x 170 x 35 mm

Basket has a Silicone mat on the inner base to secure instruments

Length Approx: 270 mm
Width Approx: 170 mm
Height Approx: 35 mm
Product Code: CBM77T



Small Instrument Basket 170 x 50 x 25 mm

Basket has a Silicone mat on the inner base to secure instruments

Length Approx: 170 mm
Width Approx: 50 mm
Height Approx: 25 mm
Product Code: CBM63T



Safety Stack System
Ultra lightweight and sturdy
Full range of sizes in stainless steel
Lids and other accessories available



Safety Stack System
Maximum Protection
Retractable Handle
Pressed Metal Rounded Corners

CBM BASKETS *"A single tray for a complete process!"*


Product Codes & Descriptions

Code	Description
CBM75/5	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 180
CBM75/3	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 90
CBM75/2	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 70
CBM75/1	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 50
CBM75	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 45
CBM75/4	PERFORATED SHEET BASKET STAINLESS STEEL 540 X 240 X 41
CBM75LID	PERFORATED LID WITH HANDLE 540 X 240
CBM74/3	PERFORATED SHEET BASKET STAINLESS STEEL 480 X 240 X 90
CBM74/2	PERFORATED SHEET BASKET STAINLESS STEEL 480 X 240 X 70
CBM74/1	PERFORATED SHEET BASKET STAINLESS STEEL 480 X 240 X 50
CBM74	PERFORATED SHEET BASKET STAINLESS STEEL 480 X 240 X 45
CBM74/4	PERFORATED SHEET BASKET STAINLESS STEEL 480 X 240 X 41
CBM74LID	PERFORATED LID WITH HANDLE 480 X 240
CBM73/3	PERFORATED SHEET BASKET STAINLESS STEEL 400 X 240 X 90
CBM73/2	PERFORATED SHEET BASKET STAINLESS STEEL 400 X 240 X 70
CBM73/1	PERFORATED SHEET BASKET STAINLESS STEEL 400 X 240 X 50
CBM73	PERFORATED SHEET BASKET STAINLESS STEEL 400 X 240 X 45

Code	Description
CBM73/4	PERFORATED SHEET BASKET STAINLESS STEEL 400 X 240 X 41
CBM73LID	PERFORATED LID WITH HANDLE 400 X 240
CBM68	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 70
CBM72/2	PERFORATED SHEET BASKET STAINLESS STEEL 300 X 240 X 70
CBM72	PERFORATED SHEET BASKET STAINLESS STEEL 300 X 240 X 50
CBM72/1	PERFORATED SHEET BASKET STAINLESS STEEL 300 X 240 X 35
CBM72LID	PERFORATED LID WITH HANDLE 300 X 240
CBM71/3	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 90
CBM71/2	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 70
CBM71/1	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 50
CBM71	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 45
CBM71/4	PERFORATED SHEET BASKET STAINLESS STEEL 240 X 240 X 41
CBM71LID	PERFORATED LID WITH HANDLE 240 X 240
CBM79/2	PERFORATED SHEET BASKET STAINLESS STEEL 650 X 240 X 70
CBM78/2	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 200
CBM78	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 150
CBM78/1	PERFORATED SHEET BASKET STAINLESS STEEL 440 X 340 X 70

INSTRUMENTS AND ENDOSCOPES SILICONE HOLDERS


- The wide range of instruments silicone holders and perforated basket organizers allows the creation of surgical procedure instruments sets in the easiest way
- The design of the instruments silicone holders permits the safe and skilful handling of the surgical instruments and endoscopes
- Thanks to its design, silicone holders can be moved rapidly from their position in order to create new configurations
- The use of silicone and stainless steel allows to withstand steam sterilisation and washing disinfection processes.
- Other shapes and dimensions can be created on request



ENDOSCOPE PERFORATED BASKETS

Endoscopes perforated baskets design represents the skilful solution for the endoscopes reprocessing.

- The body as well as the lid are manufactured entirely in perforated stainless steel AISI304 grade.
- The round corners stainless steel basket does not create any problem of scratches while using paper wrapping or pouches.
- The baskets are electrically polished. This process avoids the colour changing due to the high temperature of the sterilisation and avoids the formation of water stains.
- The new basket water draining capacity has been improved with the realisation of a smooth surface that facilitate the expulsion of the water from them
- The design of the instruments silicone holders permits the safe and skilful handling of the endoscopes
- The use of silicone and stainless steel allows to withstand steam sterilisation and washing disinfection processes.




ENDOSCOPE PERFORATED BASKETS

Endoscope Supports

Each basket includes two silicone endoscope holders that, thanks to their design, allow the use with different endoscopes dimensions.

Moreover their height allow the downward positioning of the optical fibre connection, avoiding the formation of critical water stains .



Available Dimensions

ARTICLE	LENGTH (mm)	HEIGHT (mm)	WIDTH (mm)	ØENDO (mm)
CBM66/S	450	60	70	FROM 3 TO 8
CBM66/S1	450	60	70	FROM 8 TO 15
CBM65/S	290	60	70	FROM 3 TO 8
CBM65/S1	290	60	70	FROM 8 TO 15

What type of basket do you want? Checklist!

REQUIRED ATTRIBUTES	YES	NO	NOTES
Should the weight be less than 1 kg?			
Is price or performance the most important need?			
Do you want to use the baskets in Containers?			
Do you want to use basket only in packing process?			
Do you want to use baskets only in washing process?			
Do you want to use baskets in both processes?			
DIN Sizes?			
ISO Sizes?			
European Washer Sterilisers?			
American Washers & Sterilisers?			
Woven Wire Mesh?			
Welded Wire Mesh?			
Pressed Metal?			
Combination of above?			
Mesh Size 3 mm x 3 mm?			
Mesh Size 6 mm x 6 mm?			
Very Fine Mesh?			
Reinforced Base?			
Stackable?			
Handles & Lids integral?			
Handles & Lids Retractable?			
Option for Lids?			
Polished Stainless Steel?			
Rounded Corners?			
Square Corners?			